


RUDOLF

BECAUSE WE CARE



Exquisite Suturing Experience

Laparoscopic Needle Holder

Laparoscopic Needle Holder with Ergonomic Handle



- Stable and torsion-resistant shanks
- Overload protection
- Ergonomic and light weight handle
- Axial handle – straight and curved design




Models with TC hard metal inserts



Models with needle self righting function




Luer-Lock connector for an improved cleaning of the inner tube



Fine and precise adjustment ratchet



Needle holder with ergonomic handle and overload protection, Ø 5 mm



Ordering no.

	Jaw type	Length	Curved handle	Straight handle
	Straight, TC inserts	330 mm	LP354-033	LP353-033
		450 mm	LP354-045	LP353-045
	Curved left, TC inserts	330 mm	LP354-133	LP353-133
		450 mm	LP354-145	LP353-145
	Curved right, TC inserts	330 mm	LP354-233	LP353-233
		450 mm	LP354-245	LP353-245
	Self righting, with coated jaws	330 mm	LP354-330	LP353-330
		450 mm	LP354-345	LP353-345

Ordering no.	
RZ352-001	Cap for Luer-Lock




RUDOLF Medical GmbH + Co. KG, represented worldwide, translates its belief in lasting values by fostering close customer relationships and offering innovative high technology products, hand instruments, rigid endoscopes as well as efficient services for numerous surgical disciplines.

RUDOLF

BECAUSE WE CARE

RUDOLF Medical GmbH + Co. KG
Zollerstr. 1, 78567 Fridingen, Germany
Tel. +49 7463 9956-0
sales@rudolf-med.com
www.rudolf-med.com

Distributed by:



**MITTELSTAND
GLOBAL**
HEALTH MADE
IN GERMANY