

# NASO-PHARYNGO-LARYNGOSCOPE

## OP-30

HIGH RESOLUTION IN THE SMALLEST DIAMETER

- Up to 18.000 pixels
- Detailed vision
- High resistance



producing quality

OPTOMIC flexible endoscopes provide excellent optical characteristics for the inspection of nasal passages, the pharynx and the larynx. Both fiberscopes have been designed to facilitate examination efficiency and enhance patient comfort.

## NASO-PHARYNGO-LARYNGOSCOPE **OP-30**



FIBROLUX DUO  
LIGHT SOURCE



FIBROLUX LED  
LIGHT SOURCE



FIBER OPTIC CABLE



FIBERSCOPE HOLDER



LEAKAGE TEST DEVICE



TWO-HANDED GRIP

	OP-30
Image guide system	18.000 pixels
Outer diameter	3,4mm
Flexible part diameter	3,2mm
Depth of field	2,5-50mm
Field of view	85°
Deflection angle	130°up/down
Direction of view	0°
Working length	300mm

- High definition images (up to 18.000 pixels).
- Focusing adjustment.
- ETO valve.
- Soakable in disinfectant liquid.
- Light guide connector for Storz and ACMI (optional Wolf).
- Delivered in a protective carrying case.

### OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain  
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14  
export@optomic.com • www.optomic.com



producing quality

Specifications are subject to change without notice

