

LED HEAD LAMP OP-HL15

LED TECHNOLOGY MAKES IT EASIER

- Wireless
- Comfortable
- Pure white light



OPTOMIC
producing quality

Lightweight and discreet head lamp using powerful LED technology to produce intense white light. The head lamp is wireless, which makes it very comfortable to use.

LED HEAD LAMP **OP-HL15**



LIGHT AND SLEEK DESIGN



BATTERY HOLDER



LIGHT INTENSITY ADJUSTMENT

- 3,5 hours of continuous operating time.
- Pure white light of high intensity.
- Brightness: more than 100,000 lux.
- Adjustable leather head band.
- Ion lithium rechargeable battery.
- Weight: 190 gr.
-

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
export@optomic.com • www.optomic.com



Specifications are subject to change without notice

